Agilent

Channel  
Partner

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**TABLE 1** FCM validation for CAR T cell detection

Parameters	Validation design	Calculation	Definition
Specificity	PB without CAR T cell ( $7 < S < 10$ )	/	The ability to generate a negative result for negative samples
Sensitivity LOD LLOQ	LOD: PB without CAR T cells ( $S = 15$ ); LLOQ: 293 FT cell line samples with different concentration of CAR19 positive 293FT cells (7 dilutions, performed in triplicates).	LOD: mean + 3 SD LLOQ: above the LOD, the lowest concentration where the CV is acceptable ( $CV < 30\%$ is acceptable)	The ability to correctly identify positive populations among a sample. LOD: the lowest amount of an analyte that can reliably be detected. LLOQ: the lowest amount of an analyte that can reliably be quantitated.
Accuracy <sup>a</sup>	Detect a sample with or without CAR19 T cells by both FCM and ddPCR ( $S = 188, 7^b$ )	Described by the Spearman correlation coefficient $r$ ( $r > 0.8$ is acceptable).	How close the detected value to the real true value
Precision Intra-assay Inter-assay	Intra-assay: triplicates for each sample ( $S = 7$ ); Inter-assay: three independent experiments, triplicates for each sample per experiment ( $S = 7$ ).	Intra-assay: calculated intra-CV for each sample and the precision is indicated by the mean of intra-CV and the range. Inter-assay: mark the mean for each sample in an experiment as M. For each sample, calculate the inter-CV from the three Ms. The precision is indicated by the mean of inter-CV and the range.	The reproducibility of a test Intra-assay: how close the results are when the same sample is repeated tested Inter-assay: the variation among different assay

Abbreviations: PB, peripheral blood; S, sample size.

<sup>a</sup>Due to the lack of acknowledged reference standards, true accuracy cannot be exactly determined. We intended to choose ddPCR as a “gold standard” method, and compare the consistency of the results between FCM and ddPCR. Of particular importance is that the results of ddPCR do not absolutely represent the true value.

<sup>b</sup>One hundred and eighty-eight pairs of clinical samples and seven pairs of 293FT cell line samples.

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1. Considerations for the Development of Chimeric Antigen Receptor (CAR) T Cell Products Draft Guidance for Industry. U.S. FDA. <https://www.fda.gov/media/156896/download>
2. Cheng, J, Mao, X, Chen, C, Long, X, Chen, L, Zhou, J, et al. Monitoring anti-CD19 chimeric antigen receptor T cell population by flow cytometry and its consistency with digital droplet polymerase chainreaction. *Cytometry*. 2023; 103(1): 1626. <https://doi.org/10.1002/cyto.a.24676>
3. Sarikonda, G., Mathieu, M., Natalia, M., Pahuja, A., Xue, Q., Pierog, P. L., Trampont, P. C., Decman, V., Reynolds, S., Hanafi, L. A., Sun, Y. S., Eck, S., Hedrick, M. N., Stewart, J. J., Tangri, S., Litwin, V., & Dakappagari, N. (2021). Best practices for the development, analytical validation and clinical implementation of flow cytometric methods for chimeric antigen receptor T cell analyses. *Cytometry. Part B, Clinical cytometry*, 100(1), 79–91. <https://doi.org/10.1002/cyto.b.21985>
4. Dasyam N, George P, Weinkove R. Chimeric antigen receptor T-cell therapies: Optimising the dose. *Br J Clin Pharmacol*. 2020 Sep;86(9):1678-1689. doi: 10.1111/bcp.14281. Epub 2020 Mar 24.

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