




Biosensing Instrument



MODULARITY, HIGH SENSITIVITY, FAST DETECTION



Biosensing Instrument (BI) was founded in 2004 with a vision to create highly sensitive, modular and cost effective SPR instruments for cutting edge research.



Over the past 10 years, BI has gained recognition in developing unique label-free technologies, award winning product development, unparalleled customer service, as well as funding from NIH to continue innovating new SPR technologies.

Surface Plasmon Resonance (SPR) has emerged as a powerful detection technique due to its high sensitivity and label-free capability. This optical based real-time detection method delivers superior results across diverse categories of research. BI SPR suitable for a broad range of applications such as life science, drug discovery, electrochemical analysis, food quality and safety, environmental science, and chemical sensor development.



Innovations from Biosensing Instrument

- ❖ Integrated SPR with optical microscopy
- ❖ Precision sample delivery with BI-DirectFlow™
- ❖ New measurement paradigm with analysis modules



Why SPR?

Label Free analysis

Small molecule detection

Real time binding information


High sensitivity and simple operation



View our popular
SPR Explained video

How does SPR work?


SPR is a phenomenon that occurs when polarized light hits a metal film at the interface of media with different refractive indices. Common SPR techniques excite and detect collective oscillations of free electrons (known as surface plasmons) via the Kretschmann configuration, in which light is focused onto a metal film through a glass prism and the subsequent reflection is detected. By monitoring the reflected beam shifts vs time, molecular binding events can be studied without the hassle of labels.




Biosensing Instrument SPR Microscopy

Integrated SPR with optical microscopy

Surface Plasmon Resonance Microscopy (SPRM) combines optical imaging with SPR technology, providing spatial mapping of the binding activity on living cells. A bright field light condenser illuminates the cells grown on the SPR sensor surface and it captures the bright field image of the cells. Simultaneously, the SPR light source projects its beam at its resonance angle onto the SPR sensor surface and the reflected beam is collected by the SPRM detector.




The SPRM detector records a sensorgram at each pixel of the sensing area, providing sensorgram data and SPR mapping for the entire sensing area. SPRM provides much more localized information and makes it possible to study heterogeneous surface binding and interactions of membrane proteins of single or multiple cells in their native state.




Left: bright field image of SH-EP1 human epithelial cells. Middle: SPR image of lectin proteins binding to glycoprotein receptors on SH-EP1 cells. Right: SPR sensorgrams of selected cells.


Wang et al. Nature Chemistry 2012, 4, 846-853



Affinity binding study of anti-EGFR on A431 human epidermal cells. Shown are bright field and SPR images of anti-EGFR binding to cells. KD= 371pM.




Kinetic binding study of 8 kDa Affibody peptide on A549 adenocarcinomic human alveolar basal epithelial cells. Shown are bright field and SPR images of peptides binding to cells. $k_a = 7.5 \times 10^4 \text{ M}^{-1} \text{ s}^{-1}$, $k_d = 7.2 \times 10^{-4} \text{ s}^{-1}$, $K_D = 10.2 \text{ nM}$.



Affinity binding study of WGA (Wheat-germ Agglutinin) on live CP-18821 Barrett's esophagus derived cells. Shown are bright field and SPR images of WGA binding to cells. $K_D = 70.2 \text{ nM}$.

Label-free binding activities on live cells

Quantitative mapping of binding affinity and kinetics



While the bright field image of *E. coli* O157:H7 bacteria (black dots on left image) appears to be constant, the SPR image on right fluctuates significantly. Monitoring this fluctuation, caused by the bacterial cell nanomotion, provides insight into its metabolism.

Syal et al. ACS Nano 2016, 10, 845-85




View the SPRM video



SPRm 200 Series

Light source	690 nm
Incident angles	40-76 Deg
Molecular weight cutoff	200 Da
Baseline noise	< 0.6 RU RMS (0.1 mDeg RMS)
Temperature Control Range	6°C to 50°C (10°C below ambient temperature max)
Field of view	Bright Field: 1200 x 900 um SPR: 600 x 450 um
Resolution	Bright Field & SPR: 1 um
Stage translation/rotation	3 mm x 3 mm / 360 deg
Sample injection method	Fully automated (Autosampler option) Manual
Kinetic constant	$k_a < 1 \times 10^7 \text{ M}^{-1}\text{S}^{-1}$ $k_d > 1 \times 10^{-5} \text{ S}^{-1}$
Dissociation constant	$K_D = 10^{-3} \text{ M (1 mM) to } 10^{-12} \text{ M (1 pM)}$
Outer dimension	690 (w) x 330 (h) x 340 (d) mm
Weight	23 kg
Consumables	Bare Au Cell chamber kit Autosampler
Accessories	<ul style="list-style-type: none"> Automated sample mixing and diluting Two 384, 96, 48, and 12 plate formats

Nanometer scale binding monitoring



Kinetic binding study of protein modified nanoparticles. Edge detection enhances particle detection. Time lapsed SPR images of 0.5um nanoparticles binding to modified sensor surface are useful for developing drug delivery systems.

Our customers




Our awards



Biosensing Instrument SPR Solutions

BI-4500 Series

BI-4500 offers 5 channel flow modes and delivers high quality binding response for low immobilization and small molecule (<100 Da) detection. Equipped with BI-DirectFlow™ technology, it integrates precision sample delivery with near-zero dispersion for fast kinetic measurement. Combined with precise temperature control, autosampler and many analysis modules, it provides fully automated operation.



BI-2500 Series


The BI-2500 SPR system is a compact, cost effective, 3 channel flow SPR system. It is capable of measuring small molecules (<100 Da) for the study of protein binding interactions. Its innovative modular design gives users the flexibility to choose amongst various analysis modules for biosensor development, electrochemistry SPR, and various liquid and gas phase SPR applications.



	BI-4500 Series				BI-2500 Series	
Light source	670 nm				670 nm	
No. of sample flow channels	5 channels				3 channels	
Detection speed	4 ms				4 ms	
Incident angles	40-47 Deg (gas) 67-81 Deg (liquid)				40-47 Deg (gas) 67-81 Deg (liquid)	
Molecular weight cutoff	100 Da				100 Da	
Baseline noise	< 0.06 RU RMS (0.01 mDeg RMS)				< 0.06 RU RMS (0.01 mDeg RMS)	
Temperature control range	6°C to 50°C (10°C below ambient temperature max)				N/A	
Sample injection volume	>50 µL (application dependent)				>50 µL (application dependent)	
Sample injection method	Fully automated (Autosampler option) Semi-automated				Manual	
Kinetic constant	$k_a < 1 \times 10^8 \text{ M}^{-1}\text{S}^{-1}$ $k_d > 1 \times 10^{-6} \text{ S}^{-1}$				$k_a < 1 \times 10^8 \text{ M}^{-1}\text{S}^{-1}$ $k_d > 1 \times 10^{-6} \text{ S}^{-1} *$	
Dissociation constant	$K_D = 10^{-3} \text{ M (1 mM) to } 10^{-12} \text{ M (1 pM)}$				$K_D = 10^{-3} \text{ M (1 mM) to } 10^{-12} \text{ M (1 pM)}$	
Analysis modules	BI-DirectFlow™, EC-DualFlow™, EC SPR, Gas SPR				Flow Injection, EC-DualFlow™, EC SPR, Gas SPR	
Outer dimension	355(w) x 250 (h) x 515 (d) mm				355(w) x 215 (h) x 365 (d) mm	
Weight	11.5 kg				8 kg	
Consumables	Bare Au	Divided Au	CM Dextran (Dx)	Streptavidin (SA)	Ni-NTA	
Accessories	Autosampler				Gas Diluter	
	<ul style="list-style-type: none"> Automated sample mixing and diluting for hands-free operation Holds two trays compatible with 384, 96, 48, and 12 plate formats 				<ul style="list-style-type: none"> Used with gas analysis module Built-in precision gas pump and diluter 	


BI-DirectFlow™ Technology

BI-DirectFlow™ technology delivers sample to the sensor surface with near-zero dispersion, enabling ultra-fast kinetics and high resolution binding analysis. This unique technology enables finer observation and removal of secondary effects such as bulk refractive index shift, mass limited transport, and non-specific binding. As a result, high quality data that is more representative of true molecular binding behavior can be obtained.




Life Science Applications


Kinetic Analysis of Anti-BSA/BSA




Formation of the anti-BSA/BSA immune complex is relevant to studies of the receptor site of the red blood cells. The reference subtracted binding curves for 4 channels at varying analyte concentrations with kinetic analysis fits are overlaid in red below. $k_a = 8.6 \times 10^4 \text{ M}^{-1}\text{s}^{-1}$, $k_d = 1.5 \times 10^{-4} \text{ s}^{-1}$, $KD = 1.7 \text{ nM}$.



Electrochemical SPR Applications



Redox reactions in proteins and other molecules are known to cause conformational changes in the proteins. Such conformational changes are often too small to be monitored using structural analysis techniques, but they can be studied using electrochemical SPR.




- (a) Cyclic voltammogram of cytochrome c, immobilized on a 3-mercaptopropionic acid-functionalized Au SPR sensor chip and
- (b) the simultaneously recorded SPR angle shift vs. potential.


Boussaad et al. Anal. Chem. 2000, 72, 222-226

Biomarker Expression of Carcinoma Cells

Vascular endothelial growth factor (VEGF) is a protein biomarker produced by oxygen-hungry cells to stimulate the growth of blood vessels. Certain types of tumor cells produce abnormally large amounts of VEGF to block the action of angiogenesis inhibitors. This action is termed as “angiogenic switch”, which leads to the metastasis of the tumor by providing blood supply for new tumors to grow. In this application, SKOV-3 ovarian cancer cells were cultured on the ceiling of the BI Flow Injection analysis module and VEGF secretion was measured using the SPR.




Liu et al. Anal. Chem. 2014, 86, 7305-7310



Small Molecule Applications


Carbonic anhydrase inhibitors are a class of pharmaceuticals that inhibit the activity of carbonic anhydrases (CA II). Clinically, these inhibitors have been used as antiglaucoma agents to alleviate mountain sickness and to manage neurological disorders. Acetazolamide (molecular weight = 222 Da) is one of the commonly used small molecule carbonic anhydrase inhibitors.




Binding curves (black) and simulation results (orange) of acetazolamide to pre-immobilized CAII. $KD = 14.7 \text{ nM}$, $k_a = 84.1 \times 10^6 \text{ M}^{-1}\text{s}^{-1}$, $k_d = 6.1 \times 10^{-2} \text{ s}^{-1}$.

Gas SPR Applications for Biosensor Development

BI SPR provides superior sensitivity in small molecule detection, which is critical for characterizing polymers and thin films at solid-gas interfaces. Adsorption of the small molecules onto the sensing materials, such as polymers, can be detected by observing SPR angle shift. The effect of three gases (Ethanol, Methanol and Acetone) were studied on four different sensing surfaces, resulting in quantitative characterization of the sensing polymers used for the gas sensors development.





www.BiosensingUSA.com



 +1.480.491.2777

 info2bi@BiosensingInstrument.com



進階生物科技股份有限公司

台北總公司 02-26959935
台中辦事處 04-22373160

免付費專線 0800251302
新竹辦事處 03-5640901

傳真 02-26958373
台南辦事處 06-2760235

www.level.com.tw
高屏台東辦事處 07-5539956