

Cellometer[®] K2

Image Cytometer
for Cell Counting & Analysis




PBMCs
Primary Hepatocytes
Stem Cells
Splenocytes
Tumor Suspension
and Other Primary Cells

level
Discovering
new boundaries

進階生物科技股份有限公司

服務專線：0800-251302


Nexcelom
Bioscience



Analysis of Cells from Heterogeneous Samples

- Whole Blood
- Peripheral Blood
- Cord Blood
- Bone Marrow
- Bronchoalveolar Lavage (BAL)

Primary Hepatocytes: Cell Count and Viability



PBMC Analysis in the Presence of Red Blood Cells

Measure PBMCs from whole blood without lysing. Obtain baseline PBMC concentration and viability prior to biomarker studies

Nucleated Cell Concentration & Viability

Evaluate cord blood and bone marrow samples

GFP Transfection Efficiency & Viability

Quickly and easily monitor DNA, RNA, and siRNA transfection

Analysis of Clumpy & Irregular-Shaped Cells


Nexcelom's proprietary pattern-recognition software enables accurate analysis of >98% of mammalian cell types

Cell Line Analysis

Automatically capture fluorescent cell images, concentration, Trypan blue or PI viability, and mean diameter in 60 seconds!

Cell Based Assays

- Cell Cycle
- Apoptosis
- GFP




Proven Performance in Many Research Areas



- **Clinical Immunology:** PBMCs
- **DMPK:** Primary Hepatocytes
- **Regenerative Medicine:** Stem Cells
- **Transplantation:** Nucleated Cells
- **Vaccine Development:** Splenocytes
- **Oncology:** Cell Lines, Cell Cycle, Apoptosis
- **Basic Research:** Primary Cells / Cell Lines / GFP

Optimized for
Primary Cell
Analysis




Contact Nexcelom regarding your cell type

Cellometer K2 Image Cytometer for Cell Counting & Analysis from Nexcelom Bioscience

FEATURES

- ➔ Dual FL/BR Channels
- ➔ Easily Edit and Import Assays
- ➔ Images for Data Verification
- ➔ Cell Size Histograms




- ➔ Analysis of Primary Hepatocytes
- ➔ Viability of WBC in Whole Blood
- ➔ Accurate PBMC Counts in the Presence of Red Blood Cells
- ➔ Total Nucleated Cell Count & Viability
- ➔ One-Step Cell Concentration & Viability
- ➔ Cell Cycle Analysis
- ➔ Apoptosis Analysis
- ➔ GFP Detection

ASSAYS

“ The Cellometer K2 instrument is a simple, quick, visual cell counting platform. We love the clear images and easily discerned confirmation that the cell count data is real because you can see the count for your own two eyes. It has allowed us to move away from difficult flow cytometer and subjective hemocytometer methods.


How It Works



Pipette 20 µl of Cell Sample



Insert Counting Chamber



Select Assay & Click Count

Assay: Immune cells, low RBC	
Cell Type F1: A_Immune Cells_Low RBC (AO)	
Cell Type F2: A_Immune Cells_Low RBC (PI)	
Sample ID: PBMC_AOPI_Dry demo-2	
Dilution: 2.00	
Count	Concentration
Total: 1750	6.06x10⁶ cells/mL
Live: 1662	5.75x10⁶ cells/mL
Dead: 88	3.03x10⁵ cells/mL
Mean Diameter	
7.1 micron	
7.2 microns	
5.8 micron	Viability: 95.0%

Get Results

Cellometer K2 Image Cytometer

Optimized Analysis of Primary Cells



Features of the Cellometer K2

Dual Fluorescence and Bright Field Imaging: staining of both live and dead cells in heterogeneous samples

User-Friendly Software and Assay Selection: Enhanced inter-operator reproducibility, minimal training, auto-save option

Fast Results: Obtain cell images, counts, size measurements, and viability calculations in 60 seconds

Small Sample Size: Only 20 μ l of sample

Broad Dynamic Range: Measurable concentration range of 1×10^5 to 1×10^7 cells/mL using Nexcelom's proprietary de-clustering function

Many Compatible Dyes: Trypan blue, AO, PI, EB, 7AAD, AO/PI, AO/EB, Calcein AM, CFDA, Calcein AM/PI, CFDA/PI

“ I am currently using the Cellometer K2 in our lab, mostly to count T cells, PBMCs, and tumor cells. I use them for cell culture, and later the cells are used for further assays like ELISA and FACS. The instrument is very accurate, especially with AOPI.



Advantages of Cellometer Image Cytometer

- ➔ **Cell Imaging**
 - Verify cell morphology and counted live/dead cells
 - Export cell images for presentations and publications
- ➔ **Pattern Recognition Software**
 - Accurately count cells in clumps
 - Count irregular-shaped cells
 - Eliminate debris from cell counts
 - Differentiate cells based on size
- ➔ **Automated Data Management**
 - Pre-set assays and automated reports
 - Archive sample images and auto-save results
- ➔ **Maintenance-free System**
 - Disposable counting chambers – no wash steps
 - No required instrument maintenance

Learn why thousands of users, including the top ten pharmaceutical companies, trust Cellometer.

On-Line Demonstrations are completed in just 20 to 30 minutes and provide an overview of how Cellometer works using existing images of cells that interest you.

On-Site Demonstrations are a convenient way to test a Cellometer system for a specific application. An experienced Applications Specialist will arrive at your lab for a hands-on session to test your cells and show how Cellometer can enhance your workflow.

Technical Seminars are an excellent way to introduce Cellometer systems to a lab group or collaborators in different laboratories within an organization. A trained biologist will discuss and demonstrate the capabilities and advantages of Cellometer image cytometry.

Call 978-327-5340 or E-mail info@nexcelom.com today to schedule a free demonstration or technical seminar.




Dual-Fluorescence for Primary Cell Viability in Heterogeneous Samples

Live / Dead Cell Concentration using AO / PI

Dual-Fluorescence Viability, using acridine orange (AO) and propidium iodide (PI), is the recommended method for accurate viability analysis of primary cells, such as PBMCs, splenocytes, and stem cells, in samples containing debris and unwanted non-nucleated cell types including red blood cells.


Acridine orange (AO) and propidium iodide (PI) are nuclear staining (nucleic acid binding) dyes. AO is permeable to both live and dead cells and stains all nucleated cells to generate green fluorescence. PI enters dead cells with compromised membranes and stains all dead nucleated cells to generate red fluorescence.

Because mature mammalian red blood cells do not contain nuclei, only live and dead mononuclear cells produce a fluorescent signal. There is no need to lyse red blood cells, saving time and eliminating an extra sample preparation step.




Export to FCS Express* for Flow-Like Data Output

Cell Cycle



Apoptosis



*FCS Express 4 Flow Cytometry software is a product of De Nova Software

Performance of the Cellometer K2 Image Cytometer




Figure 1. Table of results for cell concentration dynamic range

Concentration Dynamic Range

Figure 1 depicts the dynamic range for cell concentration measurements on Cellometer K2. This data set was taken on a concentration series of cultured Jurkat cell line.

Samples from 1×10^5 – 1×10^7 cells/ml can be counted without further dilution.

The %CV at each concentration was below 10%.

Viability Dynamic Range The viability dynamic range is 0 - 100% for Cellometer K2 Image Cytometer using dual fluorescence AO/PI stain.

Sample	N Value	Average Live Cell Concentration	% Viability	CV of Concentration	CV of Viability
Jurkat	24	3.61E+06	92.2%	8.9%	1.0%
Human PBMC	10	5.94E+06	96.0%	4.7%	0.5%
Mouse Splenocyte	10	1.86E+07	88.6%	5.6%	0.7%

Figure 2. Table of results for cell concentration and viability using AOPI

Consistency and Repeability The results indicate the accuracy of the Cellometer K2 instrument in assessing the viability of Jurkat cells using AOPI for cell viability. Jurkat, human PBMC, mouse splenocytes were tested at 24, 10, and 10 sample replications, respectively. The viability average was calculated and plotted. The results show the reliability and accuracy of the Cellometer K2 in measuring cell concentration and viability of mammalian cells.

“ The Cellometer K2 has drastically changed our work flow in the lab. We are able to gather cell counts in minutes rather than waiting overnight for colonies to grow on plates. It also cuts down time in the prep of plating and error in plating/counting. The amount of time that the machine has saved us is incalculable - it has allowed us to move projects along much more quickly and with confidence. - Synlogic

Cellometer Cell Counters, Cell Analysis Systems & Image Cytometry
Nexcelom offers a wide range of Cellometer systems developed and optimized for specific applications and cell types.



www.nexcelom.com/products

Cellometer®
Simply Counted



Celigo®
Image Cytometer

Which Instrument is Right for Me?

Features	Bright Field Cell Counters			Fluorescent Viability Cell Counters				Image Cytometers				
	Mini	Auto T4	Auto 1000	Auto 2000	X1	X2	K2	Vision CBA	Vision CBA (10x)	Celigo BF	Celigo 4 Channel	Celigo 5 Channel
Cell / Sample Type												
Cell Line	X	X	X	X			X	X		X	X	X
Cultured Primary Cells	X	X	X	X			X	X		X	X	X
Algae									X			
Platelets						X			X			
Low Concentration Cell Lines				X			X	X		X	X	X
Yeast (Clean Sample)					X	X			X			
Yeast (Messy Sample)						X			X			
Primary cells (Messy Sample*)				X			X	X			X	X
PBMCs, Splenocytes, Stem Cells				X			X	X			X	X
Hepatocytes							X	X			X	X
Adipocytes***				X			X	X		X	X	X
Cell-Based Assay **					X	X	X	X	X	X	X	X
Apoptosis (Annexin V-FITC/PI)							X	X	X		X	X
Apoptosis (Caspase Activity)							X	X	X		X	X
Autophagy (CytolD-green)								X	X			
Cell Proliferation (CFSE)								X	X		X	X
Cell Cycle (PI)					X	X	X	X	X		X	X
GFP Transfection				X		X	X	X	X		X	X
RFP Transfection								X	X		X	X
Mitochondrial Potential (JC-1)								X	X		X	X
Multi-drug Resistance (ABC Transporter)								X	X		X	X
Surface Marker Analysis								X	X		X	X
Vitality (Calcein-AM/PI)						X	X	X	X		X	X
Vitality (CFDA-AM)						X						
Image Cytometry**								X	X		X	X

* A messy sample is a heterogeneous sample containing unwanted cell types, such as red blood cells, in addition to the cells of interest.
** FCS Express license must be purchased in order to perform Cell Based Assay or Image Cytometry analysis
*** Cellometer CHT4-PD300 slides are required for cells greater than 80µm in diameter

Nexcelom products are for RESEARCH USE ONLY and are not approved for diagnostic or therapeutic use. © Copyright 2017 Nexcelom Bioscience LLC. All Rights Reserved. 1001235 Rev.F.11/18



進階生物科技股份有限公司

台北總公司 02-26959935
台中辦事處 04-22373160

免付費專線 0800251302
台南辦事處 06-2760235

傳真 02-26958373
高屏台東辦事處 07-5539956

新竹辦事處 03-5640901

www.level.com.tw

經銷商 榮陽長庚區:康寧 02-28200822