

THE MULTIPLEX DIGITAL PCR COMPANY



CRYSTAL DIGITAL PCR™



CRYSTAL DIGITAL PCR™

Stilla's Crystal Digital PCR™ technology for absolute quantification of nucleic acids is based on cutting-edge microfluidic technology that integrates the digital PCR process in a single consumable, reducing hands-on time and interaction, so you can simply insert your sample and walk away.

The sample is partitioned using a network of microchannels into a large array of individual droplets, also called a droplet crystals. The individual droplets are partitioned into identical sizes – each its own reaction compartment – before PCR. PCR then amplifies the fluorescent target in the droplets. The result is the accurate and precise quantification of multiple targets in a single reaction, saving you both time and precious sample.

With a combination of powerful imaging, the flexibility to use up to 6 detection channels, and proprietary software for analysis, Crystal Digital PCR™ is designed to offer an unmatched level of confidence in digital PCR measurement.

FEATURED PRODUCTS

PREPARE

Sapphire chip

Up to 12 samples / run
Up to 28,000 droplets / sample
For high sensitivity applications



Opal chip

Up to 48 samples / run
Up to 20,000 droplets / sample
For flexible throughput

AMPLIFY



Geode

Thermocycler for partitioning and amplification

naica® system > **PREPARE | AMPLIFY | READ | ANALYZE**

The naica® system harnesses the key principles of digital PCR and is designed to provide a sensitive, fast, and easy-to-use solution.



Easy-to-use digital PCR solution

On-chip integrated workflow
Minimal hands-on time



Flexible Digital PCR

Up to 6-color target multiplexing
From 12 to 48 samples / run



Fast time to results

Prepare (10 min), Amplify (2 hrs), Read (<20 min),
Analyze (5 min)


READ



Prism6

High multiplexing capacity
3- or 6-color system configuration
Wide range of fluorophore compatibility

ANALYZE



Crystal Miner software

Intuitive software for visual inspection
and analysis

CRYSTAL DIGITAL PCR™

PREPARE



Load the reaction mix into the wells of the chips.

naica® PCR MIX reagents are developed for enhanced performance for Crystal Digital PCR™ on the naica® system



AMPLIFY




Place the chips in the Geode.

Launch the Crystal Digital PCR™ program.

Thousands of droplets are created from each sample

PCR amplification is performed immediately after the Droplet Crystal is generated




READ



Transfer the chips to the Prism6.*

Droplet crystals are imaged using up to 6 fluorescent channels



* 20 minutes for two detection channels – may vary if more channels are utilized

OVERVIEW



ANALYZE





Measure the concentrations of targeted nucleic acids with Crystal Miner software.

Intuitive tool for visual inspection and analysis of Crystal Digital PCR™ experiments

Automatic quality control for experiment performance

Automatic identification of positive and negative droplets for all fluorescence channels



01
EXPLORE

Explore data
using 1D, 2D, or 3D dot
plots

02
INSPECT


Inspect droplet
crystals
and zoom in for
quality control (droplet
fluorescence versus
spatial position in the
crystal)

03
VISUALIZE

Get results
target nucleic acid
concentrations

04
ACCESS

Access raw data
and export all
experiment details and
results




CRYSTAL DIGITAL PCR™

Sapphire chip



Capacity	4 samples / chip
Number of chips per box	12
Input volume	25 μ L / sample
Number of droplets per sample	Up to 28000 droplets
Dynamic range of detection (95%)	~ 5 logs

Opal chip



Capacity	16 samples / chip
Number of chips per box	12
Input volume	7 μ L / sample
Number of droplets per sample	Up to 20,000
Dynamic range of detection (95%)	~ 5 logs

Geode



Capacity with Sapphire chip	<ul style="list-style-type: none">• Up to 12 samples (3 chips) / run• Up to 36 samples / 8h shift (3 runs)
Capacity with Opal chip	<ul style="list-style-type: none">• Up to 48 samples (3 chips) / run• Up to 144 samples / 8h shift (3 runs)
Footprint (WxDxH)	35 x 37 x 29 cm
Power supply	100-240 V~ // 50 Hz // 750W

SPECIFICATIONS



Prism6




Capacity	Up to 3 chips / run
Scan time	Less than 20 min*
Compatible fluorophores with 3-color configuration**	FAM, Atto550, Cy [®] 5...
Compatible fluorophores with 6-color configuration**	FAM, YY [®] , Atto550, ROX, Cy [®] 5, Atto700
Power Supply	110-240 V / 400W
Footprint	50 x 60x 36cm

*Dependent upon number of samples and specific dPCR assay design

**Equivalent fluorophores can be used within the same wavelengths range

Crystal Miner software



License & updates fees	Free
File format	.ncx (naica crystal experiment) .ncr (naica crystal result)
Data visualization	Droplet crystal imaging; 1D, 2D, 3D plots
File size	60 MB / well (.ncx) 10 MB / well (.ncr)
Export format	.csv / .xlsx / .png / .yaml
Threshold	Automatic and manual
Compensation matrix	Automatic and manual
Pooling	Manual

CRYSTAL DIGITAL PCR™

A RANGE OF APPLICATIONS

Experience the flexibility of the naica® system for a wide range of key nucleic acid detection and quantification applications.



To discover more
www.stillatechnologies.com



MKT-00130



進階生物科技股份有限公司

台北總公司 02-26959935 免付費專線 0800251302 傳真 02-26958373
經銷商 榮陽長庚區:康寧 02-28200822 www.level.com.tw



進階官網



FB粉絲團