# Transwell®

Permeable Supports 領導者


More than **25** Years! Over **30,000** Citations!

Create a More Natural Environment for Your Cells

1 Research Tool Used 5 Ways

## 1 Drug screening in complex cell types


Permeable supports enable the formation of a tight cell layer at the top of the membrane, which allows assessment of transport, diffusion, secretion, permeability, and drug uptake of compounds added to the cells.

**Examples**  
ADME/Tox screening (blood-brain barrier, intestinal epithelium).

## 2 Differentiation of specialized cells


Permeable supports allow cultures at the air-liquid interface: the apical side of the cells is exposed to air, while the basolateral side is immersed in liquid media. This mirrors the cells' natural environment and promotes their full differentiation.

**Examples**  
*In vitro* tissue modeling of epithelial cells (epidermis, airway epithelia, disease models, organoids).

## 5 Chemotaxis and migration assays


Permeable supports enable the analysis of the cells' ability to migrate through the membrane pores towards a chemoattractant (e.g., a tumor spheroid) grown in a Corning spheroid plate.

**Examples**  
Immune response, cancer metastasis.

## 4 Invasion assays


Permeable supports pre-coated with extracellular matrices serves as a barrier for non-invasive cells, while presenting an appropriate environment to study cell invasion.

**Examples**  
Invasion capacity of normal/malignant or compound treated/non-treated cells towards a chemoattractant.

## 3 2D, 3D, and complex co-culture studies

Permeable supports enable 2D, 3D, and co-culture of different cell types, allowing the exchange of secreted factors through membrane pores, without cell-to-cell contact.

**Examples**  
2D, 3D cell culture, intracellular communication, or cell metabolism that influences other cell types (gene expression, secretion).

5 VERSATILE USES OF TRANSWELL PERMEABLE SUPPORTS

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
實驗目的	細胞種類	建議孔徑大小 (µm)
Transport and Permeability Studies	Caco-2	0.4, 1.0
	MDCK	0.4, 1.0
Epithelial Cell Polarity	Epithelial cells	0.4
Co-culture	Stem, neuronal, and various others	0.4, 1.0
Organoid	Kidney	0.4
Tissue Engineering/ Air-Liquid Interface	Human skin model: Airway epithelial cells, disease model (e.g., COVID-19)	0.4, 3.0
Toxicity Testing	Mouse fibroblasts	3.0
	Human lung	0.4
Angiogenesis	Endothelial, HMVEC, HUVEC	3.0
Migration	Endothelial, HUVEC, HMVEC	3.0
	Neutrophils, PMNs	3.0
	Lymphocytes, macrophages, monocytes	3.0, 5.0
	Neuronal cells	3.0
	Dendritic cells	3.0, 5.0, 8.0
	Neurite outgrowth	1.0, 3.0
	Epithelial fibroblasts	8.0
Invasion	Leukocytes	3.0, 5.0
	Smooth muscle	8.0
	Melanoma	8.0
	Glioma	8.0
	Lymphoma, Jurkat	5.0, 8.0
	Osteoblasts	8.0
	Breast cancer	5.0, 8.0
Endothelial	3.0, 5.0, 8.0	

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Format	Growth Area		Recommend Medium Volume	
	Dish/Plate Well	Insert	Dish/Plate Well	Insert
100mm Dish	55 cm <sup>2</sup>	44 cm <sup>2</sup>	13 mL	9 mL
6 Well Plate	9.5 cm <sup>2</sup>	4.67 cm <sup>2</sup>	2.6 mL	1.5 mL
12 Well Plate	3.8 cm <sup>2</sup>	1.12 cm <sup>2</sup>	1.5 mL	0.5 mL
24 Well Plate	1.9 cm <sup>2</sup>	0.33 cm <sup>2</sup>	600 µL	100 µL
96 Well Plate	0.32 cm <sup>2</sup>	0.143 cm <sup>2</sup>	235 µL	75 µL

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產品功能等同 Corning Matrigel #354234

✓ E6906 - 10ML

Growth-factor reduced, without phenol red

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
\*HTS: High Throughput Screen, 自動化操作

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Cell Migration and Invasion with Transwell Video	3D Glioma Blood Brain Barrier Model for HTS of Tumoricidal Capability	Spheroid-based 3D Invasion Model for Drug Screening	Hepatocyte Differentiation from ESCs
Kidney Organoid Formation Protocol	Optimizing your Chemotaxis or Invasion Assay	Cellular Fluorescence Imaging on Transwell	Trypsinization Procedure for Transwell
			Preparation of Transwell for Histology

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