

MSC研究 鑑別抗體產品列表

Product	Cat.No. (Antibody clone ID)
---------	-----------------------------

MSC PRIMARY ANTIBODIES-POSITIVE MAKERS

CD73	11-0739-42(AD2)
CD90	11-0909-42(5E10)
CD44	11-0441-82(IM7)
CD36	50-0369-42(NL07)
CD105	12-1057-42(SN6)
Nestin	14-9843-82(10C2), 53-9843-82(10C2)
Stro-1	14-6688-82(STRO-1)

MSC PRIMARY ANTIBODIES-NEGATIVE MAKERS

CD11b	11-0113-42(CBRM1/5), 11-0118-42(ICFR44)
CD14	14-0149-82(61D3)
CD19	11-0199-42(HIB19), 17-0198-42(SJ25C1), 46-0197-42(J3-129)
CD34	11-0349-42(4H11)
CD79a	12-0792-42(HM47)
CD45	11-0459-42(HI30), 11-9459-42(2D1), 14-9457-82(CD45-2B11)
HLA-DR	11-9952-42(L243), 11-9956-42(LN3), 11-0748-42(5-329)