


ELISA 團購 揪團 一起買最划算

產品特色	Uncoated ELISA kits	Coated ELISA kits	獨家 Instant ELISA kits	推薦 ProcartaPlex immunoassays
靈敏度	<10 pg/mL	<10 pg/mL	<5 pg/mL	<10 pg/mL
Kit 偵測範圍	<5–250 pg/mL	<5–250 pg/mL	<7.8–500 pg/mL	<5–2,000 pg/mL
所需樣品量	10–100 µL	10–100 µL	25–50 µL	25–50 µL
可否偵測多重免疫因子	No	No	No	Yes
每well 測得的目標蛋白數量	1	1	1	1–80
定量表現系統	HRP-TMB (colorimetric)	HRP-TMB (colorimetric)	HRP-TMB (colorimetric)	RPE (fluorescent)
定量偵測儀器	Microplate reader	Microplate reader	Microplate reader	Luminex® instrument (Luminex® 100/200™, MAGPIX®, or FLEXMAP 3D® system)



市場獨賣 Instant ELISA kit 超強技術

- 靈敏度升級
- 樣品用量省
- 減少手作流程誤差

省下大量時間讓您輕鬆做研究!

熱賣品項 **HOT**

SKU	Product Name	Size
88-7064-88	IL-6 Mouse Uncoated ELISA Kit	10 x 96 Tests
88-7324-88	TNF alpha Mouse Uncoated ELISA Kit	10 x 96 Tests
88-7013-88	IL-1 beta Mouse Uncoated ELISA Kit	10 x 96 Tests
88-7314-88	IFN gamma Mouse Uncoated ELISA Kit	10 x 96 Tests
88-7137-88	IL-13 Mouse Uncoated ELISA Kit	10 x 96 Tests
88-7066-88	IL-6 Human Uncoated ELISA Kit	10 x 96 Tests
88-7346-88	TNF alpha Human Uncoated ELISA Kit	10 x 96 Tests
88-8086-88	IL-8 Human Uncoated ELISA Kit	10 x 96 Tests

Multiskan SkyHigh 酵素免疫分析儀

產品特色

- 6、12、24、48、96至384孔，以及各式比色管皆可做全光譜掃描的複合式光譜分析儀
- 加大觸控式螢幕
- 全波長機型200-1000 nm
- 內建DNA, RNA, protein定量模式
- 讀取迅速96孔盤6秒內完成



全面滿足
您實驗的需求



Fluoroskan FL 冷螢光分析儀

產品特色

- 擁有最齊全與高品質的濾光片選擇
激發波長範圍：320 ~ 700 nm
發散波長範圍：360 ~ 670 nm
- 支援 6~384 well 盤型，具單孔多點測量模式
- 最高可升級3支分注器，可應用於快速冷螢光反應實驗
- Top/Bottom螢光判讀模式，適合均質樣本及細胞試驗
- 高敏感度冷光讀取系統
- 產品應用：測量螢光、冷光、雙冷光、FRET、Calcium flux



Varioskan LUX 全功能雙核心微盤分析儀

產品特色

- 雙核心光學模組，全波長四光柵/高精準低背景，濾片式模組/高靈敏
- 自動抓取最佳偵測模式，不需要調整及設定PMT參數
- 內建溫控及可調式震盪功能，可升級自動分注器及氣體控制模組
- 支援多種樣品格式，6~1536 well、 μ Drop & μ Drop Duo Plates
- 激發頻寬可調：5 nm & 12 nm，為增加波長精準度
激發波長範圍：200 ~ 1000 nm
發散波長範圍：270 ~ 840 nm，可進行全波長掃描
- 產品應用：螢光、時差式螢光、吸收光、冷光、BRET、FRET、TR-FRET、AlphaScreen、細胞分析等